


Product information: 20537

impression disinfection

- suitable for all kinds of impressions, used dentures and dental workpieces
- aldehyde- and phenol-free
- with excellent cleaning effect and broad range of action
- bactericidal, fungicidal and virucidal activity (incl. TB)
- gentle to the material
- can be used in an ultrasonic bath



CE  
0481

## RHEOSEPT-ID form

Ready-to-use disinfectant for dental items with good cleaning effect

### Product description

Ready-to-use cleaning and disinfectant bath for dental workpieces and impressions made of polyethers, silicones, alginates and hydrocolloids by immersion. Alcohol-free with a pleasant scent. Range of action: bactericidal, fungicidal, virucidal / enveloped viruses (HIV, HBV, HCV, TB, vaccinia, BVDV). This medical product is CE certified (CE 0481).

### Fields of application

Disinfection of dental work pieces and impressions.

### Dosing instructions

Immerse workpiece in the product (ready to use). The intended application times are to be defined in the specific cleaning and disinfection plan. Recommended application times are shown in the application table. Recommended application times are between 5 minutes and 15 minutes, also under high load mechanics. After disinfection, rinse the

workpiece under running water. Also rinse any cavities until no more foam appears. Useful life: 14 days. If heavily soiled, early renewal of the product is recommended. To ensure the highest possible security of application it is recommended to change the used solution daily. Always comply with the recommended contact times! Product mixing must be avoided, since especially in aldehyde products precipitates, discolouration and efficacy losses are possible. Before the first use, troughs and application aids should be rinsed thoroughly with water. Not suitable for aluminium surfaces. Follow the instructions of the material manufacturer. Do not additionally heat up the bath. Suitable for ultrasonic baths.

# RHEOSEPT®

about expert reports, please contact our office staff:  
Telephone: +49 2241 - 3923-70 e-mail: info@rheosol.de.

## Environmental information

Wachendorff-Chemie GmbH produce under application of advanced, safe and environmentally friendly procedures and in compliance with high quality standards.

## Contact times

Dosing specifications	Dosing details	Remark
contact time 5 minutes	undiluted	bactericide/fungicide (acc. DGHM) high load/mechanics
contact time 30 minutes	undiluted	enveloped viruses (HBV/HIV/HCV/vacciniavirus), mycobacteria (TBC)
contact time 5 minutes	undiluted	bactericide (EN 1276), virucidal at enveloped viruses according to RKI / DVV (HIV, HBV, HCV, BVDV, vaccinia and influenza viruses)

## Composition

100g contain: 0.44g of bis(3-aminopropyl)dodecylamine, 0.14g of didecylmethyl ammonium chloride. Further ingredients: surfactants, complexing agents, corrosion protection.

## Physico-chemical data

- pH-Wert 9-10  
(Anwendung)

## Labelling in accordance with the EC-directives

- H412 - Harmful to aquatic life with long lasting effects.
- P273 - Avoid release to the environment.
- P501 - Dispose of contents/container at hazardous or special waste collection point.

## Special notes

Not subject to labelling in accordance with Directive 1999/45/EC. Note batch number and expiry date separately printed on the packaging.

## Forms of delivery

Article ID	Package size	Packaging unit
20537352	RHEOSEPT-ID form bottle 1000 ml	carton
20537050	RHEOSEPT-ID form canister 5 litres	canister

## Registrations

This medical product is CE certified (CE 0481).

## Certificate and information

For additional information on this product, please see → [www.rheosol.com](http://www.rheosol.com) – For specific questions or information

**RHEOSOL®**