


Product information: 20546

drill bath

- suitable for rotary instruments
- aldehyde-free
- with excellent cleaning effect and broad range of action
- bactericidal, fungicidal and virucidal activity (incl. TB)
- gentle to the material
- can be used in an ultrasonic bath



CE
0481

RHEOSEPT-ID drill

Ready-to-use drill bath with good cleaning effect

Product description

Ready-to-use cleaning and disinfecting bath for steel and diamond drills, root canal instruments and other stainless steel instruments. Other areas of application in the field of foot care, such as nail nippers/clippers, scissors, drills, cutting and grinding instruments. High corrosion protection, material compatibility and cleaning power. Pleasant scent. Range of action: bactericidal, fungicidal, virucidal / enveloped viruses (HIV, HBV, HCV, TB, vaccinia, BVDV). This medical product is CE certified (CE 0481).

Fields of application

For medical instruments, particularly in the field of dentistry and foot care. Particularly suitable for rotary instruments such as dental drills, milling and grinding equipment.

Dosing instructions

Immerse instruments in the product (ready to use). The

intended application times are to be defined in the specific cleaning and disinfection plan. Recommended application times are shown in the application table. Recommended application times are between 5 minutes and 15 minutes, also under high load mechanics. After disinfection, rinse the instrument under running water. Also rinse any cavities until no more foam appears. Useful life: 14 days. If heavily soiled, early renewal of the product is recommended. To ensure the highest possible security of application it is recommended to change the used solution daily. Always comply with the recommended contact times! Product mixing must be avoided, since especially in aldehyde products precipitates, discolouration and efficacy losses are possible. Before the first use, troughs and application aids should be rinsed thoroughly with water. Not suitable for aluminium surfaces. Follow the instructions of the material manufacturer. Do not additionally heat up the bath. Suitable for ultrasonic baths.

about expert reports, please contact our office staff:
Telephone: +49 2241 - 3923-70 e-mail: info@rheosol.de.

Environmental information

Wachendorff-Chemie GmbH produce under application of advanced, safe and environmentally friendly procedures and in compliance with high quality standards.

Contact times

Dosing specifications	Dosing details	Remark
contact time 5 minutes	undiluted	bactericide/fungicide (acc. DGHM) high load/mechanics
contact time 30 minutes	undiluted	enveloped viruses (HBV/HIV/HCV/vacciniavirus), mycobacteria (TBC)
contact time 5 minutes	undiluted	bactericide (EN 1276), virucidal at enveloped viruses according to RKI / DVV (HIV, HBV, HCV, BVDV, vaccinia and influenza viruses)

Composition

100g contain: 0.44g of bis(3-aminopropyl)dodecylamine, 0.14g of didecylmethyl ammonium chloride. Further ingredients: surfactants, complexing agents, corrosion protection.

Physico-chemical data

- pH-value 9-10 (application)

Labelling in accordance with the EC-directives

- H412 - Harmful to aquatic life with long lasting effects.
- P273 - Avoid release to the environment.
- P501 - Dispose of contents/container at hazardous or special waste collection point.

Special notes

Not subject to labelling in accordance with Directive 1999/45/EC. Note batch number and expiry date separately printed on the packaging.

Forms of delivery

REF	Package size	Packaging unit
20546352	RHEOSEPT-ID drill bottle 1000 ml	carton
20546050	RHEOSEPT-ID dent Bohrerbad Kanister 5 Liter	canister

Registrations

This medical product is CE certified (CE 0481).

Certificate and information

For additional information on this product, please see → www.rheosol.com – For specific questions or information