

Product information: 20502

surface disinfection

- with excellent cleaning effect and broad range of action
- protects and cares
- highly economical due to low application concentrations
- aldehyde-free
- excellent bactericidal, fungicidal and viricidal activity with optimal application capabilities
- effective against influenza, noro and MRSA
- surface disinfectant as high concentrate
- VAH listed



CE
0481

RHEOSEPT-FD plus

Liquid surface disinfection concentrate without the addition of aldehyde

Product description

Disinfectant for the institutional area on the basis of cationic active ingredients, without the addition of aldehyde, for surface disinfection, wiping disinfection as well as disinfecting surface cleaning. Listed with the VAH (Verband für Angewandte Hygiene e.V.) acc. to the guidelines of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.). Adequate for water-resistant surfaces. Spectrum of activity: bactericidal, fungicidal, viricidal (HIV, HBV, HCV, noro, vaccinia, BVDV). This medical product is CE certified (CE 0481).

Fields of application

For application in hospitals, kindergartens, schools, foster homes, swimming pools and other public facilities for disinfection.

Dosing instructions

For disinfection prepare the dilution (0,5%-60 Min., 1%-15 Min., 1%-30 Min., 3%-5 min.) depending on the intended residence time and apply it thinly and wiping on the pre-cleaned surface. According to the application area of the medical product the surface is to rinse or wipe carefully after disinfection in order to remove product residues. Assure complete wetness, the disinfectant is still effective after long residence times. Always use fresh prepared dilutions.

Spectrum of activity

Spectrum of activity: bactericidal, fungicidal, viricidal (HIV, HBV, HCV, noro, vaccinia, BVDV).

Contact times

Dosing specifications	Dosing details	Remark
contact time 60 minutes	0,5	bactericide/fungicide (VAH) high load/mechanics
contact time 30 minutes	1%	bactericide/fungicide (VAH) high load/mechanics
contact time 15 minutes	1%	bactericide/fungicide (VAH) high load/mechanics
contact time 5 minutes	3 %	bactericide/fungicide (VAH) high load/mechanics
contact time 60 minutes	1%	enveloped viruses (HIV, HBV, HCV, BVDV, vaccinia, influenza) EN 14476
contact time 5 minutes	3 %	enveloped viruses (HIV, HBV, HCV, BVDV, vaccinia, influenza) EN 14476
contact time 60 minutes	1%	norovirus (MNV) EN 14476

Composition

100g contain: 5.0g didecyldimethylammoniumchloride; 1.9g quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl-, chlorides; 4.0g N-(3-Aminopropyl)-N-dodecylpropane-1.3-diamine. Further ingredients: complexing agents, alcohol, scents, linalool, hexyl cinnamal, limonenes.

Physico-chemical data

- pH-value (application) approx. 9
- pH-value (concentrate) 11.0-11.5

Labelling in accordance with the EC-directives

- H302 - Harmful if swallowed.
- H314 - Causes severe skin burns and eye damage.
- H411 - Toxic to aquatic life with long lasting effects.
- P260 - Do not breathe dust/fume/gas/mist/vapours/spray.
- P280 - Wear protective gloves/protective clothing/eye protection/face protection.
- P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P314 - Get medical advice/attention if you feel unwell.
- P101 - If medical advice is needed, have product container or label at hand.

Specific advices

According to the application area of the medical product the surface is to rinse or wipe carefully after disinfection in order to remove product residues. Note batch number and expiry date separately printed on the packaging. Expiry date: 36 months as of manufacturing date.

Forms of delivery

REF	Package size	Packaging unit
20502352	RHEOSEPT-FD plus bottle 1000 ml	carton
20502050	RHEOSEPT-FD plus canister 5 litres	canister
20502100	RHEOSEPT-FD plus canister 10 litres	canister

Accessories

REF	Package size	Packaging unit
48029002	Bottle holder for 4 bottles à 1000 ml	piece
48020002	Canister pump plug 51mm, 30 ml pump stroke canister pump for 5 and 10-litres-canister blue or white with UN-code Riser pipe length 300 mm	piece
48017002	Canister tap for 5 and 10-litres-canister with plug 51	piece
48029003	Canister hook for 10 litres canister for wall fastening	piece
48041001	Protective glasses with integrated side protection and plastic glasses in acc. to EN 166 single-packed in bags, carton with 100 pieces	piece
48038002	Canister key 51+61 mm	ID7966
48036001	Protective gloves yellow, size S, package with 12 pairs 100 % natural latex, lined with cotton, rolled edging, suitable for the food processing industry	package
48036006	Protective gloves yellow, size M, package with 12 pairs 100 % natural latex, lined cotton, rolled edging, suitable for the food processing industry	package
48036005	Protective gloves yellow, size L, package with 12 pairs 100 % natural latex, lined cotton, rolled edging, suitable for the food processing industry	package
48036003	Protective gloves yellow, size XL, package with 12 pairs 100 % natural latex, lined cotton, rolled edging, suitable for the food processing industry	package
48036002	Protective gloves blue, 1 pair Latex, size 10, acid and alkaline-resistant (30 cm long) package with 12 pairs	package

Certificate and information

For additional information on this product, please see → www.rheosol.com – For specific questions or information about expert reports, please contact our office staff: Telephone: +49 2241 - 3923-70 e-mail: info@rheosol.de.

Admission

Listed with the VAH (Verbund für Angewandte Hygiene e.V.) acc. to the guidelines of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.). This medical product is

CE certified (CE 0481).

Environmental information

Wachendorff-Chemie GmbH produce under application of advanced, safe and environmentally friendly procedures and in compliance with high quality standards.