

Product information: 24603

hand decontamination

- skin-friendly and nurturing washing product with germ-reducing effect
- tested acc. to DIN EN 1499 - Hygienic handwash
- suitable washing product in the medical field
- skin neutral ph-value and without colourants
- Fragrance-free formulation for the food processing industry.
- with skin-nurturing allantoin
- skin-friendly - was awarded highest marks when dermatologically tested



RHEOSEPT-HD wash

Hygienic washing lotion for hand disinfection

Product description

Hygienic washing lotion for quick and easy hand and skin hygiene. Skin-friendly and nurturing washing product with germ-reducing effect. Tested acc. to EN 1499. Good decontaminating effect even when diluted with water. Contains allantoin for the care of stressed skin.

Fields of application

The product is suitable for hygienic handwashing in the food industry and other sensitive areas.

Dosing instructions

Use 3 ml of undiluted product on moistened hand. Rub in well for 30 seconds and rinse under running water. The application can be used several times a day. We recommend to dry hands with disposable towels.

Contact times

Dosing specifications	Dosing details	Remark
contact time 30 seconds	3ml undiluted	hygienic hand washing EN 1499

Composition

In 100 g contain: 0.3 g polyhexamethylene biguanide (PHMB). Other excipients: Allantoin.

Physico-chemical data

- pH-value (application) 6.0

Labelling in accordance with the EC-directives

- H319 - Causes serious eye irritation.
- P301 + P330 + P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
- P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P337 + P313 - If eye irritation persists: Get medical advice/attention.
- EUH208 - Contains polyhexamethylene biguanide hydrochloride. May produce an allergic reaction.

Certificate and information

For additional information on this product, please see → www.rheosol.com – For specific questions or information about expert reports, please contact our office staff: Telephone: +49 2241 - 3923-70 e-mail: info@rheosol.de.

Admission

This medical product is CE-certified.

Environmental information

Wachendorff-Chemie GmbH produce under application of advanced, safe and environmentally friendly procedures and in compliance with high quality standards.

Forms of delivery

Article ID	Package size	Packaging unit
24603353	RHEOSEPT-HD wash system cartridge 500 ml	carton
24603355	RHEOSEPT-HD wash system cartridge 1000 ml	carton
24603060	RHEOSEPT-HD wash canister 5 litres	carton
24603100	RHEOSEPT-HD wash canister 10 litres	canister

Accessories

Article ID	Description	Packaging unit
48017001	canister tap plug 45 canister tap	piece
48017006	Canister tap plug 40 Canister tap for 5 litres-lightweight canister	piece
48026032	Euro-dosing pumps as cap for Euro system cartridges Bag with 10 dosing pumps Dispenser cap for Euro bottle when using mobile without dispenser Dosing pump with special sealing and additional riser pipe	bag
48020001	Canister pump plug 45 mm 30 ml pump stroke canister pump for standard canister/top speed-canister 10 litres	piece
48038004	Canister key plug 45 Canister	piece
48029003	Canister hook for 10 litres canister for wall fastening	piece